

The image shows a video player interface. The main content area is a video frame with a black border. The video content features a blue, crinkled paper background. In the center, there is a small, partially visible image of a person's hands packing items into a bag. Overlaid on this background is the title 'Packing a Prep Kit' in a large, red, serif font with a white outline. At the bottom right of the video frame, there are two small, teal-colored buttons: one labeled 'Menu' and another with a left-pointing chevron symbol.

Packing a Prep Kit

Menu





Items required

Menu



3 Medium Sized Gloves



Menu



4 x 10cm x 10cm Gauze



Menu



A Towel



Menu



Alice Forceps



Menu



Mouth Mirror



Menu



150cc Metal Cup



Menu



Autoclave Bags 13.25" x 20"

This provides a Large Sterile Surface to work from when opened

Autoclave Bags 10.5" x 16"



Menu



Autoclave Bags 10.5" x 16"



Menu



Autoclave Bags 10.5" x 16"

This Smaller Bag, with the Prepping Kit inside, fits into the Autoclave more easily

Collect together the following Items

- Medium Size Gloves
- 10cm x 10cm / 4'' x 4''
Gauze Squares
- A Towel

Menu



Collect together the following Items

- 150cc Metal Bowl
- A Pair of Alice Forceps
- A Mouth Mirror

Menu



Collect together the following Items

- Autoclave Bags 13.25" x 20"
- Autoclave Bags 10.5" x 16"
- Betadyne Solution

Menu



Fold the Towel into a Small Square



Menu



Fold the Towel into a Small Square



Menu



Fold the Towel into a Small Square



Menu



Take Four 10cm x 10cm
Gauze Squares



Menu



**Put One on Top of Each Other
on the Worksurface**



Menu



**Place the Folded Towel
on the Gauze**



Menu



**Place the Metal Cup
on the Towel**



Menu



Take 3 Medium Sized Gloves



Menu



**Fold the Cuffs over back
on each other**



Menu



**Put the Folded Gloves
on Top of Each Other**



Menu



Put the Gloves on Top of the Cup



Menu



Put a Toothbrush on the Right of the Gloves



Menu



Alice Forceps and Gauze



Menu



Put an Alice Forceps on the Right of the Toothbrush



Menu



Put a Mouth Mirror on the Right of the Alice Forceps



Menu



Take a Large Autoclave Bags
13.25" x 20"



Menu



Place the Prepping Items inside



Menu



**Fold the End of the Bag over
along the Taped Surface**



Menu



**Fold the End of the Bag over
along the Taped Surface**



Menu



Tear off the Tape



Menu



**Press the Sticky Side of the
Taped Area down firmly**



Menu



The Sealed First Bag



Menu



**Fold the Edges of the Bag
back on Themselves**



Menu



**Fold the Edges of the Bag
back on Themselves**



Menu



**Folding the Bag this way
makes the Contents Visible**



Menu



The Bag ready to be placed in the Smaller Bag



Menu



Autoclave Bags 10.5" x 16"



Menu



**Place the Bag
inside the Smaller Bag**



Menu



**Place the Bag
inside the Smaller Bag**



Menu



The Bag ready to seal



Menu



**Fold the End of the Bag over
along the Taped End**



Menu





Seal the Bag



Menu



Seal the Bag



Menu



Seal the Bag



Menu



Seal the Bag



Menu



The Sealed Bag



Menu



The Completed Prep Kit



Menu





The End

Menu

